

Recombinant Human ATRIP

Catalog No: #GP12614



Package Size: #GP12614-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ATRIP
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 592-791 amino acids of human ATRIP
Target Name	ATR interacting protein
Accession No.	Swissprot:Q8WXE1Gene Accession:BC030597
Uniprot	Q8WXE1
GeneID	84126;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes an essential component of the DNA damage checkpoint. The encoded protein binds to single-stranded DNA coated with replication protein A. The protein also interacts with the ataxia telangiectasia and Rad3 related protein kinase, resulting in its accumulation at intranuclear foci induced by DNA damage. Multiple transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only