

Recombinant Human BRD7

Catalog No: #GP12639



Package Size: #GP12639-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BRD7
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 452-651 amino acids of human BRD7
Target Name	bromodomain containing 7
Other Names	BP75; NAG4; CELTIX1
Accession No.	Swissprot:Q9NPI1Gene Accession:BC050728
Uniprot	Q9NPI1
GeneID	29117;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein which is a member of the bromodomain-containing protein family. The product of this gene has been identified as a component of one form of the SWI/SNF chromatin remodeling complex, and as a protein which interacts with p53 and is required for p53-dependent oncogene-induced senescence which prevents tumor growth. Pseudogenes have been described on chromosomes 2, 3, 6, 13 and 14. Alternative splicing results in multiple transcript variants.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only