

Recombinant Human BCAN

Catalog No: #GP12642



Package Size: #GP12642-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BCAN
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 712-911 amino acids of human BCAN
Target Name	brevican
Other Names	BEHAB; CSPG7
Accession No.	Swissprot:Q96GW7Gene Accession:BC022938
Uniprot	Q96GW7
GeneID	63827;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the lectican family of chondroitin sulfate proteoglycans that is specifically expressed in the central nervous system. This protein is developmentally regulated and may function in the formation of the brain extracellular matrix. This protein is highly expressed in gliomas and may promote the growth and cell motility of brain tumor cells. Alternate splicing results in multiple transcript variants.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only