

Recombinant Human BRF2

Catalog No: #GP12643



Package Size: #GP12643-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BRF2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 220-419 amino acids of human BRF2
Target Name	BRF2, RNA polymerase III transcription initiation factor 50 kDa subunit
Other Names	BRFU; TFIIIB50
Accession No.	Swissprot:Q9HAW0 Gene Accession:BC010648
Uniprot	Q9HAW0
GeneID	55290;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes one of the multiple subunits of the RNA polymerase III transcription factor complex required for transcription of genes with promoter elements upstream of the initiation site. The product of this gene, a TFIIIB-like factor, is directly recruited to the TATA-box of polymerase III small nuclear RNA gene promoters through its interaction with the TATA-binding protein.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only