

Recombinant Human BRF1

Catalog No: #GP12644



Package Size: #GP12644-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BRF1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	BRF1, RNA polymerase III transcription initiation factor 90 kDa subunit
Other Names	BRF; CFDS; hBRF; BRF-1; GTF3B; TAF3C; TAF3B2; TF3B90; TAFIII90; TFIIB90; HEL-S-76p
Accession No.	Swissprot:Q92994Gene Accession:BC016743
Uniprot	Q92994
GeneID	2972;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes one of the three subunits of the RNA polymerase III transcription factor complex. This complex plays a central role in transcription initiation by RNA polymerase III on genes encoding tRNA, 5S rRNA, and other small structural RNAs. The gene product belongs to the TF2B family. Several alternatively spliced variants encoding different isoforms, that function at different promoters transcribed by RNA polymerase III, have been identified.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only