

Recombinant Human BRPF3

Catalog No: #GP12646



Package Size: #GP12646-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BRPF3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1006-1205 amino acids of human BRPF3
Target Name	bromodomain and PHD finger containing 3
Accession No.	Swissprot:Q9ULD4Gene Accession:BC117387
Uniprot	Q9ULD4
GeneID	27154;
Storage	-20~-80°C, pH 7.6 PBS

Background

Component of the MOZ/MORF complex which has a histone H3 acetyltransferase activity.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only