

# Recombinant Human BTF3L4

Catalog No: #GP12650



Package Size: #GP12650-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human BTF3L4
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	basic transcription factor 3-like 4
Accession No.	Swissprot:Q96K17Gene Accession:BC022371
Uniprot	Q96K17
GeneID	91408;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The function of this protein remains unknown.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only