

## Recombinant Human LRIF1

Catalog No: #GP12667



Package Size: #GP12667-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human LRIF1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 570-769 amino acids of human LRIF1
Target Name	ligand dependent nuclear receptor interacting factor 1
Other Names	RIF1; C1orf103
Accession No.	Swissprot:Q5T3J3Gene Accession:BC008115
Uniprot	Q5T3J3
GeneID	55791;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Represses the ligand-induced transcriptional activity of retinoic acid receptor alpha (RARA). This repression may occur through direct recruitment of histone deacetylases.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only