

TADA2A Antibody

Catalog No: #43114

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TADA2A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total TADA2A protein.
Immunogen Description	Synthetic peptide of human TADA2A
Target Name	TADA2A
Other Names	ADA2; ADA2A; KL04P; hADA2; TADA2L
Accession No.	Swiss-Prot#: O75478Gene ID: 6871
Uniprot	O75478
GeneID	6871;
Concentration	1.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Gastric cancer tissue using #43114 at dilution 1/100.



Immunohistochemical analysis of paraffin-embedded Human Esophagus cancer tissue using #43114 at dilution 1/100.

Background

Many DNA-binding transcriptional activator proteins enhance the initiation rate of RNA polymerase II-mediated gene transcription by interacting functionally with the general transcription machinery bound at the basal promoter. Adaptor proteins are usually required for this activation, possibly to acetylate and destabilize nucleosomes, thereby relieving chromatin constraints at the promoter. The protein encoded by this gene is a transcriptional activator adaptor and has been found to be part of the PCAF histone acetylase complex. Several alternatively spliced transcript variants encoding different isoforms of this gene have been described, but the full-length nature of some of these variants has not been determined.?

Note: This product is for in vitro research use only