

Recombinant Human C1orf112

Catalog No: #GP12671



Package Size: #GP12671-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human C1orf112
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 654-853 amino acids of human C1orf112
Target Name	chromosome 1 open reading frame 112
Accession No.	Swissprot:Q9NSG2Gene Accession:BC107168
Uniprot	Q9NSG2
GeneID	55732;
Storage	-20~-80°C, pH 7.6 PBS

Background

The function of this protein remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only