

Recombinant Human TSEN15

Catalog No: #GP12674



Package Size: #GP12674-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TSEN15
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	tRNA splicing endonuclease subunit 15
Other Names	sen15; C1orf19
Accession No.	Swissprot:Q8WW01Gene Accession:BC022030
Uniprot	Q8WW01
GeneID	116461;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a subunit of the tRNA splicing endonuclease, which catalyzes the removal of introns from tRNA precursors. Alternative splicing results in multiple transcript variants. There is a pseudogene of this gene on chromosome 17.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only