

Recombinant Human C1R

Catalog No: #GP12676



Package Size: #GP12676-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human C1R
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 506-705 amino acids of human C1R
Target Name	complement C1r subcomponent
Accession No.	Swissprot:P00736Gene Accession:BC035220
Uniprot	P00736
GeneID	715;
Storage	-20~-80°C, pH 7.6 PBS

Background

C1r B chain is a serine protease that combines with C1q and C1s to form C1, the first component of the classical pathway of the complement system.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only