

# Recombinant Human GATD3B

Catalog No: #GP12679



Package Size: #GP12679-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human GATD3B
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 42-268 amino acids of human GATD3B
Target Name	glutamine amidotransferase like class 1 domain containing 3B
Other Names	KNP-I
Accession No.	Swissprot:A0A0B4J2D5Gene Accession:BC003587
Uniprot	A0A0B4J2D5
GeneID	8209;
Storage	-20~-80°C, pH 7.6 PBS

## Background

May play a role in promoting mitochondrial enlargement in cone photoreceptor cells.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only