

Recombinant Human C21orf59

Catalog No: #GP12680



Package Size: #GP12680-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human C21orf59
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	chromosome 21 open reading frame 59
Other Names	Kur; FBB18; CILD26; C21orf48
Accession No.	Swissprot:P57076Gene Accession:BC000709
Uniprot	P57076
GeneID	56683;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene a protein plays a critical role in dynein arm assembly and motile cilia function. Mutations in this gene result in primary ciliary dyskinesia.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only