

Recombinant Human C3orf38

Catalog No: #GP12684



Package Size: #GP12684-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human C3orf38
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 130-329 amino acids of human C3orf38
Target Name	chromosome 3 open reading frame 38
Accession No.	Swissprot:Q5JPI3Gene Accession:BC024188
Uniprot	Q5JPI3
GeneID	285237;
Storage	-20~-80°C, pH 7.6 PBS

Background

May be involved in apoptosis regulation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only