

Recombinant Human NADK2

Catalog No: #GP12686



Package Size: #GP12686-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NADK2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 164-442 amino acids of human NADK2
Target Name	NAD kinase 2, mitochondrial
Other Names	DECRD; MNADK; NADKD1; C5orf33
Accession No.	Swissprot:Q4G0N4Gene Accession:BC062567
Uniprot	Q4G0N4
GeneID	133686;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a mitochondrial kinase that catalyzes the phosphorylation of NAD to yield NADP. Mutations in this gene result in 2,4-dienoyl-CoA reductase deficiency. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only