

ADAMTS13 Antibody

Catalog No: #43116

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ADAMTS13 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Immunogen Description	Synthetic peptide of human ADAMTS13
Target Name	ADAMTS13
Other Names	VWFCP; C9orf8; vWF-CP; ADAM-TS13; ADAMTS-13
Accession No.	Swiss-Prot#: Q76LX8Gene ID: 11093
Uniprot	Q76LX8
GeneID	11093;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #43116 at dilution 1/20.



Immunohistochemical analysis of paraffin-embedded Human colorectal cancer tissue using #43116 at dilution 1/20.

Background

This gene encodes a member of a family of proteins containing several distinct regions, including a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin type 1 (TS) motif. The enzyme encoded by this gene specifically cleaves von Willebrand Factor (vWF). Defects in this gene are associated with thrombotic thrombocytopenic purpura. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only