

Recombinant Human CA6

Catalog No: #GP12690



Package Size: #GP12690-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CA6
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 18-179 amino acids of human CA6
Target Name	carbonic anhydrase 6
Other Names	CA-VI; GUSTIN
Accession No.	Swissprot:P23280Gene Accession:BC034350
Uniprot	P23280
GeneID	765;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is one of several isozymes of carbonic anhydrase. This protein is found only in salivary glands and saliva and protein may play a role in the reversible hydration of carbon dioxide though its function in saliva is unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only