

Recombinant Human CABP5

Catalog No: #GP12692



Package Size: #GP12692-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CABP5
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	calcium binding protein 5
Other Names	CABP3
Accession No.	Swissprot:Q9NP86Gene Accession:BC126133
Uniprot	Q9NP86
GeneID	56344;
Storage	-20~-80°C, pH 7.6 PBS

Background

The product of this gene belongs to a subfamily of calcium binding proteins, which share similarity to calmodulin. Calcium binding proteins are an important component of calcium mediated cellular signal transduction. Expression of this gene is retina-specific. The mouse homolog of this protein has been shown to express in the inner nuclear layer of the retina, suggested its role in neuronal functioning. The specific function of this gene is unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only