

Recombinant Human CALCOCO2

Catalog No: #GP12700



Package Size: #GP12700-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CALCOCO2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 247-446 amino acids of human CALCOCO2
Target Name	calcium binding and coiled-coil domain 2
Other Names	NDP52
Accession No.	Swissprot:Q13137Gene Accession:BC015893
Uniprot	Q13137
GeneID	10241;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a coiled-coil domain-containing protein. The encoded protein functions as a receptor for ubiquitin-coated bacteria and plays an important role in innate immunity by mediating macroautophagy. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only