

Recombinant Human CALHM1

Catalog No: #GP12701



Package Size: #GP12701-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CALHM1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 147-346 amino acids of human CALHM1
Target Name	calcium homeostasis modulator 1
Other Names	FAM26C
Accession No.	Swissprot:Q8IU99Gene Accession:BC036193
Uniprot	Q8IU99
GeneID	255022;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a calcium channel that plays a role in processing of amyloid-beta precursor protein. A polymorphism at this locus has been reported to be associated with susceptibility to late-onset Alzheimer's disease in some populations, but the pathogenicity of this polymorphism is unclear.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only