

Recombinant Human CANT1

Catalog No: #GP12705



Package Size: #GP12705-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CANT1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 202-401 amino acids of human CANT1
Target Name	calcium activated nucleotidase 1
Other Names	DBQD; DBQD1; SCAN1; SHAPY; SCAN-1
Accession No.	Swissprot:Q8WVQ1Gene Accession:BC017655
Uniprot	Q8WVQ1
GeneID	124583;
Storage	-20~-80°C, pH 7.6 PBS

Background

This protein encoded by this gene belongs to the apyrase family. It functions as a calcium-dependent nucleotidase with a preference for UDP. Mutations in this gene are associated with Desbuquois dysplasia with hand anomalies. Alternatively spliced transcript variants have been noted for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only