

Recombinant Human CPM

Catalog No: #GP12708



Package Size: #GP12708-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CPM
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 244-443 amino acids of human CPM
Target Name	carboxypeptidase M
Accession No.	Swissprot:P14384Gene Accession:BC022276
Uniprot	P14384
GeneID	1368;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a membrane-bound arginine/lysine carboxypeptidase. Its expression is associated with monocyte to macrophage differentiation. This encoded protein contains hydrophobic regions at the amino and carboxy termini and has 6 potential asparagine-linked glycosylation sites. The active site residues of carboxypeptidases A and B are conserved in this protein. Three alternatively spliced transcript variants encoding the same protein have been described for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only