

Recombinant Human CCDC137

Catalog No: #GP12718



Package Size: #GP12718-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CCDC137
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	coiled-coil domain containing 137
Other Names	RaRF
Accession No.	Swissprot:Q6PK04Gene Accession:BC009369
Uniprot	Q6PK04
GeneID	339230;
Storage	-20~-80°C, pH 7.6 PBS

Background

The function of this protein remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only