

Recombinant Human CCDC113

Catalog No: #GP12721



Package Size: #GP12721-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CCDC113
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 178-377 amino acids of human CCDC113
Target Name	coiled-coil domain containing 113
Other Names	HSPC065
Accession No.	Swissprot:Q9H0I3Gene Accession:BC126445
Uniprot	Q9H0I3
GeneID	29070;
Storage	-20~-80°C, pH 7.6 PBS

Background

Component of centriolar satellites contributing to primary cilium formation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only