

# Recombinant Human HMMR

Catalog No: #GP12728



Package Size: #GP12728-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human HMMR
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 395-594 amino acids of human HMMR
Target Name	hyaluronan mediated motility receptor
Other Names	CD168; IHABP; RHAMM
Accession No.	Swissprot:O75330Gene Accession:BC108905
Uniprot	O75330
GeneID	3161;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is involved in cell motility. It is expressed in breast tissue and together with other proteins, it forms a complex with BRCA1 and BRCA2, thus is potentially associated with higher risk of breast cancer. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only