

## Recombinant Human SLC3A2

Catalog No: #GP12734



Package Size: #GP12734-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SLC3A2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 431-630 amino acids of human SLC3A2
Target Name	solute carrier family 3 member 2
Other Names	4F2; CD98; MDU1; 4F2HC; 4T2HC; NACAE; CD98HC
Accession No.	Swissprot:P08195Gene Accession:BC001061
Uniprot	P08195
GeneID	6520;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene is a member of the solute carrier family and encodes a cell surface, transmembrane protein. The protein exists as the heavy chain of a heterodimer, covalently bound through di-sulfide bonds to one of several possible light chains. The encoded transporter plays a role in regulation of intracellular calcium levels and transports L-type amino acids. Alternatively spliced transcript variants, encoding different isoforms, have been characterized.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only