

Recombinant Human CDH19

Catalog No: #GP12742



Package Size: #GP12742-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CDH19
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 246-445 amino acids of human CDH19
Target Name	cadherin 19
Other Names	CDH7; CDH7L2
Accession No.	Swissprot:Q9H159Gene Accession:BC015877
Uniprot	Q9H159
GeneID	28513;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is one of three related type II cadherin genes situated in a cluster on chromosome 18. The encoded protein is a calcium dependent cell-cell adhesion glycoprotein containing five extracellular cadherin repeats. Loss of cadherins may be associated with cancer formation. Alternative splicing results in multiple transcript variants for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only