

Recombinant Human CDK2AP2

Catalog No: #GP12743



Package Size: #GP12743-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CDK2AP2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	cyclin dependent kinase 2 associated protein 2
Other Names	p14; DOC-1R
Accession No.	Swissprot:O75956Gene Accession:BC002850
Uniprot	O75956
GeneID	10263;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that interacts with cyclin-dependent kinase 2 associated protein 1. Pseudogenes associated with this gene are located on chromosomes 7 and 9. Alternatively spliced transcript variants have been observed for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only