

Recombinant Human CDKAL1

Catalog No: #GP12745



Package Size: #GP12745-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Recombinant Human CDKAL1 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Full length fusion protein |
| Target Name | CDK5 regulatory subunit associated protein 1 like 1 |
| Accession No. | Swissprot:Q5VV42Gene Accession:BC064145 |
| Uniprot | Q5VV42 |
| GeneID | 54901; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

The protein encoded by this gene is a member of the methylthiotransferase family. The function of this gene is not known. Genome-wide association studies have linked single nucleotide polymorphisms in an intron of this gene with susceptibility to type 2 diabetes.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only