

Recombinant Human CEP70

Catalog No: #GP12756



Package Size: #GP12756-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CEP70
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-212 amino acids of human CEP70
Target Name	centrosomal protein 70
Other Names	BITE
Accession No.	Swissprot:Q8NHQ1Gene Accession:BC016050
Uniprot	Q8NHQ1
GeneID	80321;
Storage	-20~-80°C, pH 7.6 PBS

Background

Plays a role in the organization of both preexisting and nascent microtubules in interphase cells. During mitosis, required for the organization and orientation of the mitotic spindle.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only