

## Recombinant Human DPF3

Catalog No: #GP12760



Package Size: #GP12760-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human DPF3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 158-357 amino acids of human DPF3
Target Name	double PHD fingers 3
Other Names	CERD4; BAF45C
Accession No.	Swissprot:Q92784Gene Accession:BC060801
Uniprot	Q92784
GeneID	8110;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the D4 protein family. The encoded protein is a transcription regulator that binds acetylated histones and is a component of the BAF chromatin remodeling complex. Alternate splicing results in multiple transcript variants encoding different isoforms.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only