

# Recombinant Human CCT8L2

Catalog No: #GP12761



Package Size: #GP12761-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CCT8L2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 358-557 amino acids of human CCT8L2
Target Name	chaperonin containing TCP1 subunit 8 like 2
Other Names	CESK1
Accession No.	Swissprot:Q96SF2Gene Accession:BC033797
Uniprot	Q96SF2
GeneID	150160;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Possible molecular chaperone; assists the folding of proteins upon ATP hydrolysis.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only