

# Recombinant Human CHI3L2

Catalog No: #GP12770



Package Size: #GP12770-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CHI3L2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 191-390 amino acids of human CHI3L2
Target Name	chitinase 3 like 2
Other Names	CHIL2; YKL39; YKL-39
Accession No.	Swissprot:Q15782Gene Accession:BC011460
Uniprot	Q15782
GeneID	1117;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is similar to bacterial chitinases but lacks chitinase activity. The encoded protein is secreted and is involved in cartilage biogenesis. Several transcript variants encoding different isoforms have been found for this gene.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only