

Recombinant Human CHMP6

Catalog No: #GP12775



Package Size: #GP12775-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CHMP6
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-201 amino acids of human CHMP6
Target Name	charged multivesicular body protein 6
Other Names	VPS20
Accession No.	Swissprot:Q96FZ7Gene Accession:BC010108
Uniprot	Q96FZ7
GeneID	79643;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the chromatin-modifying protein/charged multivesicular body protein family. Proteins in this family are part of the ESCRT-III (endosomal sorting complex required for transport III) which degrades surface receptors, and in biosynthesis of endosomes.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only