

Recombinant Human CIAO1

Catalog No: #GP12780



Package Size: #GP12780-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CIAO1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 140-339 amino acids of human CIAO1
Target Name	cytosolic iron-sulfur assembly component 1
Other Names	CIA1; WDR39
Accession No.	Swissprot:O76071Gene Accession:BC001395
Uniprot	O76071
GeneID	9391;
Storage	-20~-80°C, pH 7.6 PBS

Background

Key component of the cytosolic iron-sulfur protein assembly (CIA) complex, a multiprotein complex that mediates the incorporation of iron-sulfur cluster into extramitochondrial Fe/S proteins. Seems to specifically modulate the transactivation activity of WT1. As part of the mitotic spindle-associated MMXD complex it may play a role in chromosome segregation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only