

Recombinant Human CLASP1

Catalog No: #GP12788



Package Size: #GP12788-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CLASP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1390-1538 amino acids of human CLASP1
Target Name	cytoplasmic linker associated protein 1
Other Names	MAST1
Accession No.	Swissprot:Q7Z460Gene Accession:BC032563
Uniprot	Q7Z460
GeneID	23332;
Storage	-20~-80°C, pH 7.6 PBS

Background

CLASPs, such as CLASP1, are nonmotor microtubule-associated proteins that interact with CLIPs (e.g., CLIP170; MIM 179838). CLASP1 is involved in the regulation of microtubule dynamics at the kinetochore and throughout the spindle (Maiato et al., 2003 [PubMed 12837247]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only