

# Recombinant Human CLUL1

Catalog No: #GP12793



Package Size: #GP12793-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CLUL1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 267-466 amino acids of human CLUL1
Target Name	clusterin like 1
Other Names	RA337M
Accession No.	Swissprot:Q15846Gene Accession:BC025381
Uniprot	Q15846
GeneID	27098;
Storage	-20~-80°C, pH 7.6 PBS

## Background

CLUL1 is a retina-specific clusterin-like protein which shows nearly 25% identity to clusterin at the protein level. It is expressed specifically in cone photoreceptor cells.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only