

# Recombinant Human CMPK1

Catalog No: #GP12794



Package Size: #GP12794-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CMPK1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	cytidine/uridine monophosphate kinase 1
Other Names	CK; CMK; UMK; CMPK; UMPK; UMP-CMPK
Accession No.	Swissprot:P30085Gene Accession:BC014961
Uniprot	P30085
GeneID	51727;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes one of the enzymes required for cellular nucleic acid biosynthesis. This enzyme catalyzes the transfer of a phosphate group from ATP to CMP, UMP, or dCMP, to form the corresponding diphosphate nucleotide. Alternate splicing results in both coding and non-coding transcript variants.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only