

Recombinant Human CNDP2

Catalog No: #GP12795



Package Size: #GP12795-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CNDP2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 276-475 amino acids of human CNDP2
Target Name	CNDP dipeptidase 2 (metallopeptidase M20 family)
Other Names	CN2; CPGL; PEPA; HsT2298; HEL-S-13
Accession No.	Swissprot:Q96KP4Gene Accession:BC001375
Uniprot	Q96KP4
GeneID	55748;
Storage	-20~-80°C, pH 7.6 PBS

Background

CNDP2, also known as tissue carnosinase and peptidase A (EC 3.4.13.18), is a nonspecific dipeptidase rather than a selective carnosinase (Teufel et al., 2003 [PubMed 12473676]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only