

Recombinant Human CNOT2

Catalog No: #GP12798



Package Size: #GP12798-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CNOT2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 341-540 amino acids of human CNOT2
Target Name	CCR4-NOT transcription complex subunit 2
Other Names	NOT2; CDC36; NOT2H; HSPC131
Accession No.	Swissprot:Q9NZN8Gene Accession:BC002597
Uniprot	Q9NZN8
GeneID	4848;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a subunit of the multi-component CCR4-NOT complex. The CCR4-NOT complex regulates mRNA synthesis and degradation and is also thought to be involved in mRNA splicing, transport and localization. The encoded protein interacts with histone deacetylases and functions as a repressor of polymerase II transcription. Alternatively spliced transcript variants have been observed for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only