

Recombinant Human COG4

Catalog No: #GP12801



Package Size: #GP12801-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human COG4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 586-785 amino acids of human COG4
Target Name	component of oligomeric golgi complex 4
Other Names	COD1; CDG2J
Accession No.	Swissprot:Q9H9E3Gene Accession:BC013347
Uniprot	Q9H9E3
GeneID	25839;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a component of an oligomeric protein complex involved in the structure and function of the Golgi apparatus. Defects in this gene may be a cause of congenital disorder of glycosylation type IIj. Two transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only