

# Recombinant Human COMMD1

Catalog No: #GP12806



Package Size: #GP12806-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human COMMD1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-190 amino acids of human COMMD1
Target Name	copper metabolism domain containing 1
Other Names	MURR1; C2orf5
Accession No.	Swissprot:Q8N668Gene Accession:BC022046
Uniprot	Q8N668
GeneID	150684;
Storage	-20~-80°C, pH 7.6 PBS

## Background

COMMD1 is a regulator of copper homeostasis, sodium uptake, and NF-kappa-B (see MIM 164011) signaling (de Bie et al., 2005 [PubMed 16267171]).

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only