

Recombinant Human COPS7A

Catalog No: #GP12812



Package Size: #GP12812-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human COPS7A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-275 amino acids of human COPS7A
Target Name	COP9 signalosome subunit 7A
Other Names	CSN7; CSN7A; SGN7a
Accession No.	Swissprot:Q9UBW8Gene Accession:BC011789
Uniprot	Q9UBW8
GeneID	50813;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a component of the COP9 signalosome, an evolutionarily conserved multi-subunit protease that regulates the activity of the ubiquitin conjugation pathway. Alternatively spliced transcript variants that encode the same protein have been described.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only