

# Recombinant Human COPZ1

Catalog No: #GP12813



Package Size: #GP12813-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human COPZ1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein(177FULL)
Target Name	coatamer protein complex subunit zeta 1
Other Names	COPZ; CGI-120; HSPC181; zeta-COP; zeta1-COP
Accession No.	Swissprot:P61923Gene Accession:BC002849
Uniprot	P61923
GeneID	22818;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a subunit of the cytoplasmic coatamer protein complex, which is involved in autophagy and intracellular protein trafficking. The coatamer protein complex is comprised of seven subunits and functions as the coat protein of coat protein complex (COP)I-vesicles. Alternative splicing results in multiple transcript variants.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only