

Recombinant Human CORO7

Catalog No: #GP12818



Package Size: #GP12818-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CORO7
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 726-925 amino acids of human CORO7
Target Name	coronin 7
Other Names	CRN7; POD1; 0610011B16Rik
Accession No.	Swissprot:P57737Gene Accession:BC117289
Uniprot	P57737
GeneID	79585;
Storage	-20~-80°C, pH 7.6 PBS

Background

Coronins, such as CORO7, constitute an evolutionarily conserved family of WD-repeat actin-binding proteins. CORO7 plays a role in Golgi complex morphology and function.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only