

Recombinant Human CPNE4

Catalog No: #GP12827



Package Size: #GP12827-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CPNE4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 358-557 amino acids of human CPNE4
Target Name	copine 4
Other Names	CPN4; COPN4
Accession No.	Swissprot:Q96A23Gene Accession:BC014396
Uniprot	Q96A23
GeneID	131034;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene belongs to the highly conserved copine family. It encodes a calcium-dependent, phospholipid-binding protein, which may be involved in membrane trafficking, mitogenesis and development. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only