

Recombinant Human CPNE7

Catalog No: #GP12829



Package Size: #GP12829-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CPNE7
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 434-633 amino acids of human CPNE7
Target Name	copine 7
Accession No.	Swissprot:Q9UBL6Gene Accession:BC035334
Uniprot	Q9UBL6
GeneID	27132;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the copine family, which is composed of calcium-dependent membrane-binding proteins. The gene product contains two N-terminal C2 domains and one von Willebrand factor A domain. The encoded protein may be involved in membrane trafficking. Two alternatively spliced transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only