

Recombinant Human CRABP1

Catalog No: #GP12834



Package Size: #GP12834-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CRABP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-137 amino acids of human CRABP1
Target Name	cellular retinoic acid binding protein 1
Other Names	RBP5; CRABP; CRABP1; CRABP-I
Accession No.	Swissprot:P29762Gene Accession:BC022069
Uniprot	P29762
GeneID	1381;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a specific binding protein for a vitamin A family member and is thought to play an important role in retinoic acid-mediated differentiation and proliferation processes. It is structurally similar to the cellular retinol-binding proteins, but binds only retinoic acid at specific sites within the nucleus, which may contribute to vitamin A-directed differentiation in epithelial tissue.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only