

Recombinant Human CRB3

Catalog No: #GP12836



Package Size: #GP12836-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CRB3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 27-123AA amino acids of human CRB3
Target Name	crumbs 3, cell polarity complex component
Accession No.	Swissprot:Q9BUF7 Gene Accession:BC018409
Uniprot	Q9BUF7
GeneID	92359;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the Crumbs family of proteins. This gene is widely expressed in epithelial tissues where the encoded protein isoforms play various roles such as the control of cytokinesis and ciliogenesis or the formation of tight junctions. Alternative splicing results in multiple transcript variants encoding different isoforms.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only